



RAZOR[®] Sim10 Training Assay Card

Yeast and Bg Target Assay

10

Not for Diagnostic Use

Kit Part No: PATH-ASY-0091

⚠ Handle all unknown samples and waste materials as if they were capable of transmitting infectious agents.

This reagent kit is designed to simulate The 10[®] Target Kit by testing one unknown liquid or powder sample against two training assays.

Organism	Common Name	Organism	Common Name
<i>Bacillus globigii</i> / <i>atrophaeus</i>	Bg	<i>Saccharomyces cerevisiae</i>	Yeast

Kit Contents

5 mL Reagent Grade Water (1)	5 mL Unknown Sample Bottles (1)	3 mL Syringes with Cannula Tips (2)	Transfer Pipettes Pack (1)
Sample Swab Pack (1)	Instruction Booklet (1)	Freeze-dried Reagents in a Foil Bag (1)	

Run Protocol Barcode

- Verify the Run Protocol is loaded onto the instrument before preparing the pouch.
- The Run Protocol may be loaded by scanning the square protocol barcode to the right into the instrument.

NOTE: *Times, temperatures, and results analyses are included in the protocol barcode.*



Kit Part Number: PATH-ASY-0091
Protocol Code: SIMUL10A

10



RAZOR® Sim10 Training Assay Card

Yeast and Bg Target Assay

10

Pouch Barcode

NOTE: Prepare sample before loading.

- a. After the square protocol barcode has been scanned, scan the rectangular barcode on the pouch fitment. If the rectangular barcode is damaged, use the generic one shown here.



SIMUL10A-SIMUL002

NOTE: The barcode may only be used for one run. To re-use it delete the previous run.


Sample Loading Order

1 Inlet Port 1

Load 2.0 mL **Unknown Sample** with 3 mL syringe

2 Inlet Port 2

Load 0.4 mL **Reagent Grade Water** from bottle with 3 mL syringe



Inserting the Pouch into the RAZOR

- a. Insert the prepared pouch into the instrument with the label facing up and plungers toward the front of the instrument.

NOTE: The pouch fitment should fit perfectly in the groove over the insert slot.

