



BioFire® COVID-19 Test v1.1 Quick Guide

For use with FilmArray® 2.0 and FilmArray® Torch Systems



DO NOT DISCARD: Important product-specific information


For in vitro diagnostic use under Emergency Use Authorization (EUA) only

LOWER RESPIRATORY


REF 423745 (6 pk)
423744 (30 pk)

IVD Rx Only

Package Contents




Materials required but not provided



NOTE: Use clean gloves and other Personal Protective Equipment (PPE) when performing this procedure.


Step 1: Prepare Pouch

- Insert pouch into Pouch Loading Station.
- Place **Sample Injection Vial** into **red well**.
- Place **Hydration Injection Vial** into **blue well**.



Step 2: Hydrate Pouch


- Unscrew **Hydration Injection Vial** from cover.
- Remove **Hydration Injection Vial**, leaving **blue plastic cover** in Pouch Loading Station.
- Insert **Hydration Injection Vial** into **pouch hydration port**.
- Push down to puncture seal and wait as **Hydration Solution** is drawn into the pouch.



NOTE: Verify the pouch has been hydrated.

Step 3: Prepare Lower Respiratory Sample Mix

- Hold Sample Buffer tube tip facing up and firmly pinch at textured plastic tab on side of tube until seal snaps.
NOTE: Do not touch tube tip.
- Dispense Sample Buffer into **Sample Injection Vial** using a slow, forceful squeeze followed by a second squeeze.
NOTE: Avoid generating excessive foam.
- Use the Sample Swab to stir the entire specimen for **~10** seconds.
- Place the swab end into the **Sample Injection Vial** then break off at the scored breakpoint. Discard the swab handle into an appropriate waste container and close **Sample Injection Vial** lid tightly.
- Invert the Sample Injection Vial **3** times then return to **red well** of Pouch Loading Station.



WARNING: Sample Buffer is harmful if swallowed and can cause serious eye damage and/or skin irritation.



Step 4: Load Lower Respiratory Sample Mix


- Unscrew **Sample Injection Vial** from **red plastic cover**.
- Wait for **5 seconds**, then lift **Sample Injection Vial**, leaving **red plastic cover** in Pouch Loading Station.

NOTE: Waiting **5 seconds** decreases the contamination risk.

- Insert **Sample Injection Vial** into **pouch sample port**.
- Push down to puncture seal and wait as Sample Mix is drawn into the pouch.

NOTE: Verify the sample has been loaded.

LOWER RESPIRATORY



Step 5: Run Pouch

- Discard the **Sample Injection Vial** and the **Hydration Injection Vial** in a biohazard sharps container.
- Follow instructions on computer for starting a test.

Step 6: Review Results

Run Summary - Displays information about the sample and a summary of the Internal Controls and test results.

1 Internal Controls:

- If 'Passed', results are valid.
- If 'Failed' or 'Invalid', **RETEST SAMPLE** and refer to *Instructions for Use*.

Result Summary - Displays overall SARS-CoV-2 test result on the first line, followed by each individual assay result.

- ### 2
- If overall 'SARS-CoV-2' test result is 'Detected' or 'Not Detected', report the results.
 - If 'Invalid', **RETEST SAMPLE** and refer to *Instructions for Use*.


Run Details - Displays information about the pouch, protocol, run status, operator, pouch serial number, instrument, and pouch lot number.

3 Run Status:

- If 'Completed', run is complete.
- If 'Incomplete', 'Aborted', or any other error message, **RETEST SAMPLE** and refer to *Instructions for Use*.

NOTE: Refer to *Instructions for Use* for reporting information.

If repeated error messages are obtained, contact *BioFire Defense Technical Support*.



BioFire®

COVID-19 Test v1.1

BIO FIRE

www.BioFireDefense.com

Run Summary

Sample ID: Example Report

Detected: None

Run Date: 31 Dec 2019 8:00 AM

Internal Controls: Passed

Result Summary

Viruses

Not DetectedSARS-CoV-2

Not DetectedSARS-CoV-2a

Not DetectedSARS-CoV-2c

Not DetectedSARS-CoV-2d

Not DetectedSARS-CoV-2e

Not DetectedSARS-CoV-2f

Not DetectedSARS-CoV-2g

Not DetectedSARS-CoV-2h

Run Details

Pouch: COVID-19 Test v1.1

Run Status: Completed

Serial No.: 01234567

Lot No.: 012345

Protocol: Sample v3.2

Operator: Anonymous

Instrument: FA0000

Conditions of Authorization

The BioFire COVID-19 Test has not been FDA cleared or approved but has been authorized for emergency use by FDA under an EUA for use by authorized laboratories.

The BioFire COVID-19 Test has been authorized only for the detection of nucleic acid from SARS-CoV-2, not for any other viruses or pathogens.

The emergency use of the BioFire COVID-19 Test is only authorized for the duration of the declaration that circumstances exist justifying the authorization of emergency use of in vitro diagnostics for detection and/or diagnosis of COVID-19 under Section 564(b)(1) of the Federal Food, Drug and Cosmetic Act, 21 U.S.C. § 360bbb-3(b)(1), unless the declaration is terminated or authorization is revoked sooner.