




For In Vitro Diagnostic Use.

NOTE: For use with BioFire Global Fever Panel (DFA2-ASY-0004) and should not be used with any other kit.


REF DFA2-ASY-0006 IVD
CONTROL **R_x Only**



To avoid contamination, always wear clean gloves and work behind a protective shield.

Step 1: Prepare Pouch


- Insert pouch into Pouch Loading Station.
- Place **Sample Injection Vial** into **red well**.
- Place **Hydration Injection Vial** into **blue well**.



Step 2: Hydrate Pouch

- Unscrew **Hydration Injection Vial** from cover.
- Remove **Hydration Injection Vial**, leaving **blue cover** in Pouch Loading Station.
- Insert **Hydration Injection Vial** into **hydration port**.
- Push down to puncture seal and wait as **Hydration Solution** is drawn into the pouch.


NOTE: Verify the pouch has been hydrated.



Step 3: Prepare External Control

Add Sample Buffer (NOTE: There are 2 possible designs of the Sample Buffer Ampoule)

- Hold the Sample Buffer Ampoule with the tip facing up.
NOTE: Do not touch the tip of the ampoule.
- Firmly pinch textured plastic tab on the side of the ampoule until the seal snaps **or** if there is no textured tab on the side, gently twist off the plastic tab on the tip.
- Dispense Sample Buffer into Control Injection Vial using a slow, forceful squeeze, followed by a 2nd squeeze.
NOTE: Avoid generating excessive foam.
- Tightly close lid and invert the Control Injection Vial 3 times.
- Return Control Injection Vial to **red well** of Pouch Loading Station.




WARNING: Sample Buffer is harmful if swallowed and can cause serious eye damage and/or skin irritation.





Step 4: Load External Control

- Unscrew Control Injection Vial from **red cover**.
 - Wait for **5 seconds**, then remove Control Injection Vial, leaving **red cover** in Pouch Loading Station.
- NOTE:** *Waiting 5 seconds decreases the contamination risk.*
- Insert Control Injection Vial into **pouch sample port**.
 - Push down to puncture seal, then wait as control material is drawn into the pouch.




Step 5: Run Pouch

- Screw vials back into covers in Pouch Loading Station before disposing of them in a biohazard container.
 - Remove pouch from Pouch Loading Station and load into the instrument.
 - Follow instructions on screen for starting a test.
- NOTE:** *Select either **Positive External Control** or **Negative External Control** protocol.*

Step 6: Review Report


Run Information - Displays information about the sample, protocol, pouch, instrument, run status, results, and recommended action.

- Internal Controls:**
 - If 'Passed', results are valid.
 - If 'Failed' or 'Invalid', **RETEST EXTERNAL CONTROL** once.
- Run Status:**
 - If 'Completed', run is complete.
 - If 'Incomplete', 'Aborted', or any other error message, **RETEST EXTERNAL CONTROL** once.
- External Control Results:**
 - If 'Passed', report the results.
 - If 'Failed':
 - Positive External Control: **RETEST EXTERNAL CONTROL** once.
 - Negative External Control: Decontaminate the area and **RETEST EXTERNAL CONTROL** once.
 - If 'Invalid', **RETEST EXTERNAL CONTROL** once.



BIOFIRE® SHIELD™

GF Panel - Positive External Control



www.BioFireDefense.com

Run Information

Sample ID	test pouch	Run Date	12 Mar 2021 12:00 AM
Protocol	Positive External Control v3.1	Serial No.	01234567
Pouch Type	GF Panel - IVD v2.1	Lot No.	012345
Internal Controls	Passed	Operator	Anonymous
Run Status	Completed	Instrument	FA0000

1

3

Passed

Report the Results.

NOTE: Refer to the *BIOFIRE SHIELD Control Kit for the BioFire Global Fever Panel Instructions for Use* for panel specific run reports and additional information. If repeated error messages are obtained, contact *BioFire Defense Technical Support*.