




For Research Use Only. Not for use in diagnostic procedures.

Information Sheet




Viruses

Chikungunya virus
 Crimean-Congo hemorrhagic fever virus
 Dengue virus
 Ebola virus
 Lassa virus
 Marburg virus
 West Nile virus
 Yellow fever virus
 Zika virus



Bacteria

Bacillus anthracis
Francisella tularensis
Leptospira spp.
Salmonella enterica serovar Typhi
Salmonella enterica serovar Paratyphi A
Yersinia pestis



Protozoa

Leishmania spp.
Plasmodium spp.
P. falciparum
P. vivax/ovale



FilmArray Global Fever Panel - RUO with BioFire FilmArray® 2.0 System

1
Test

19 Targets. All in Under an Hour.

The RUO FilmArray Global Fever Panel tests for a comprehensive list of 9 viral, 6 bacterial and 4 protozoan targets directly from whole blood. With just one test you can identify these targets in under an hour with only 2 minutes of hands-on time.

- **Fully Automated:** Sample prep, amplification, identification, and reporting.
- **Single Instrument Integration:** Minimal equipment and consumables.
- **Freeze-dried Reagents:** Room temperature stable.

Kit Specifications

Kit	FilmArray Global Fever Panel - RUO
Part No:	DFA2-ASY-0003
Sample Types	Whole Blood
Sample Volume	~0.2 mL (200 µL)
Tests per Kit	6 Pack Kit
Reagents	
<ul style="list-style-type: none"> • Freeze-dried in durable plastic pouches • Room temperature storage 	

The purchase of these products includes a limited, nontransferable instrument license under specific claims of one or more U.S. patents as listed on BioFire Defense's website (<http://biofiredefense.com/legal/notice>) (the "Website") and owned by the University of Utah Research Foundation and/or BioFire. Any kits (i) sold with these products and/or discussed herein may be covered by one or more of the U.S. patents, as listed on the Website for the product and (ii) sold herein include a limited, nontransferable license to use the enclosed amount(s) in such kits according to the specified protocols.